## Walk-through

This walk-through is intended to be an easy to follow step-by-step guide for a "standard clinical electron" set-up, design and order electron apertures. It should be noted that the values provided are not intended to represent clinical guidelines and in no way be considered a recommendation for values that are appropriate for clinical use. The intent is that after completing this walk-through the user will have gained the knowledge an understanding of many of the details and much of the information to complete a clinical electron set up utilizing the decimal3D App. With this information the user will be able to to complete the process using the clinical protocols and requirements for your facility to set-up a clinical electron treatment, design the electron aperture and then order the electron aperture from .decimal.

- 1. Have the physician outline the treatment area in the exam room
- 2. Open the decimal3D app on the iPad and login with user name and password



- 3. Bring the patient to the simulator or treatment room
- 4. Select New Patient in the top right hand of the Patient Management Screen

		PATIENT MANAGEMENT			+ NEW PATIENT		
®	Patient Management	Q Search					
		FULL NAME	MRN	DOB		STATUS	ACTIONS
⊚ ⊖		Jane	1345	10/06/1978	506320	Planning	OPEN
		Jane doe	9876	10/10/2009	B11489	Planning	OPEN
		Keltt cup	56789		X07281	Ordered	OPEN
		Kyle	98764	05/15/1950	Q12161	Planning	OPEN
		Linda	789	09/15/2019	P10851	Planning	OPEN
		Marty Arm	6554378	07/04/2005	008003	Planning	OPEN
		Mary arm	456788	07/07/1996	109605	Ordered	OPEN
		Mary foot	6789	09/16/2019	ciralg	Planning	OPEN
	decimal 💴	October	7890	10/08/1998	010158	Planning	OPEN
	Martin D						

decimal App Documentation - https://apps.dotdecimal.com/

- 5. Enter Patient Name
- 6. Enter Date of Birth from drop dwon
- 7. Select patient Gender from drop down
- 8. Enter MRN (Medical Record Number)
- 9. The user may select the enter the Physician if desired
- 10. Select OK that the patient is entered successfully
- 11. In the top right hand corner select New Scan
- 12. Select the New Scan button on the right hand side
- 13. Place the area to be scanned in the box. The area to be scanned should appear in all red. This will require the user to stand approximately a half meter away from the patient. A scan be obtained by either

×

×

×

×

×

×

×

×

From: https://apps.dotdecimal.com/ - decimal App Documentation

Permanent link: https://apps.dotdecimal.com/doku.php?id=decimal3d:instructions\_for\_use:customer\_walkthrough&rev=1571927684

Last update: 2021/07/29 18:21