Organization Configuration

The ElectronRT App allows for users with certain permissions to view and edit organization level configuration settings. The Organization Configuration block is located in the main page of the app and is only visible to users with Physics or higher level permissions.

	▼
Patier	it Search
Seard	n for patients in the system.
► Impor	t Patients
Select	DICOM patients to import.
▼ Organ	ization Configuration
▼ Orga	nization Settings
Organiza	tion Level Settings
Add or ea	lit settings for your entire organization, including:
- CT Ov - Electr	erride materials on Device Materials
Import	t
▶ Ехро	rt Logs
Site Facili	ty Settings
Select a S	ite ID to edit settings for a facility
Site ID:	+

Fig. 1: Organization Configuration UI

Organization Settings

The Organization Settings block allows the user to view and edit settings that affect the entire organization. These settings include organization name, PDF report logo, CT override materials, and electron device materials.



Fig. 2: Organization Settings

While editing the settings in the Organization Configuration block, the UI is changed such that the user cannot click on other sections of the UI until confirming edits with "Done" or canceling edits with "Cancel". Some settings are grouped into tabs, such as the "General" and "Materials" tabs in the Organization Settings.

		Dura Ortan		
	Materials	Dose Colors		
Organ	ization Info			
	Name:	.decimal LLC		
Report Cu	stomization			
	Report Logo:		Browse	
		recommended size: 220 x 50px		

Fig. 3: Organization Edit UI

Materials

The Materials definitions allow specifying the device materials and CT override materials. The materials are represented as json in the following example format:

```
Device Material
```

```
{
    "bolus": {
        "density": 0.920,
        "name": "BlueWax",
        "relative_scattering_power": 0.920,
        "relative_stopping_power": 0.920
}
```

}

Override Material Example

```
{
    "Water": {
        "density": 1,
        "name": "Water",
        "relative_scattering_power": 1.0,
        "relative_stopping_power": 1.0
    }
}
```

Export Logs

The Export Logs block allows the user to export a file containing logs of user activity within the app. This .csv file contains data exported from the app database that keeps track of critical user activity, including (but not limited to): opening of patients and plans, plan approvals, report/DICOM exports, and hardware ordering.



Fig. 6: Export Logs UI

Site Facility Settings

The Site Facility Settings section of the Organization Configuration allows users to view and edit settings for each site in the organization. After selecting a Site ID, the site level settings are divided into five categories: Site Settings, DICOM Settings, CT Curves, QA Options, and Machine Settings.



Site Settings

The Site Settings block allows for the viewing and editing of miscellaneous site level settings including the site address, physicians, and treatment sites. In order to edit physicians, the user must first click on the physician name on the list of physicians on left side UI and then click on "Edit Physician" on the right side UI (this is also the case with treatment sites). Users can add or remove physicians and treatment sites using the left side UI.

Recipient Factored	- 0 ×
decinal an	Admin lingum Rows (1960) 1240 😝 🛽
Organization Configuration	Sin Letings
Organization Settings	Paulights Barto 108
+ Experilage	Bigang Lidens
tion famility larges	1 (voljuten elem)
Salari silari il la nali anlingo la silanily Salari 12145 (decimal) -	Physician Ser
* Six Sellings	D ann
Rysteine Altre strong dustant	Trainert Sie
	Name Hansenhalt Andrean Name Handelsen Stadt das Andrean Stadt (na Andrean
	and the second s
Testimet Elles Aller ment heterides	
THE REAL PROPERTY AND INCOME.	
0mi	
Posta a	
+ DCDB belogs	
+ CFGarees	
+ CR Cyllone	
+ Marilan Latings	



DICOM Settings

The DICOM Settings block allows the user to view and edit settings related to the export of DICOM files. These settings include the default DICOM export directory and a list of DICOM server AE titles.



Fig. 10: DICOM Settings UI

CT Curves

The CT Curves block allows for the viewing and editing of CT conversion curves. These curves must contain data for both relative stopping and scattering power. The relative stopping power data for the selected curve is displayed on a graph on the right side UI.



Fig. 11: CT Curves UI

The CT Curves are defined as json as shown in the below example (note: this data is provided as a user guide reference only and should not be used in a commissioned treatment planning system; as such, it does not contain full values of CT curve data):

```
CT Curve Example
```

```
"value": 0.0010
      },
      {
          "key": -706.0,
          "value": 0.2920
      },
   ],
   "stopping_power_curve": [
      {
          "key": -10000.0,
         "value": 0.0010
      },
      {
          "key": -1000.0,
         "value": 0.0010
      },
   ]
}
```

QA Options

The QA Options block allows the user to view and edit settings related to Dose QA. These settings include the dimensions of the water phantom, dose grid spacing, and the default measurement plane depth.

Second Extra Electroniti			- 6	з ж
Denitor Ven				
Jectrial an				
Paden Search Search Search Hoper Failures Search Sections Search Sections Search SECMULARIAN & Hight	GA Options Content GA Plantes Date: Teacort (C) Impact spotted Cater Sector States Setter Part Cater Sector (C) Care Cater Sector (C) Care	94 43.50 96 10		
Organization Configuration	THE OWNER AND A DESCRIPTION OF	12.42.43.00		
Organization Settings	Webs Pharton Stat X			
+ Experilage	Web Places Status	500 5		
Naci - State -				
* Six Settings				
+ DEDB Sellings				
+ CT Games				
OR Options Ont Options Ont the College In the sameticals				
+ Multier Letings				

Fig. 13: QA Options UI

Machine Settings

The Machine Settings block allows the user to view and edit machines for the selected site. The machine settings are divided into five tabs: General, Geometry, Applicators, Commissioning, and Advanced.

General Machine Settings

The General tab of the Machine Settings block allows for the viewing and/or editing of general machine information such as the machine name, description, serial, type, and physical SAD.

animai Restored		- D K
Spin Yas		
.decinal an		itheric Veryne Vinera 🕮 🖬 🖬
	Description (Section Company) Advance	
Padent Samh Samh in te spiers	No for Viscolin Notes Serve Se	
★ Report Patients Salest Intillitypatients to report	Bergies Australia. Sett. 1997	
* Organization Configuration	The state	
+ Dynamics Integr	LacModectly Allebringer Reve LacIModecDate 2020418131710	
 Expertisp 		
Den Facility Serlings Reset a time for a structure for a lasting		
Sill 1010 placement +		
+ Ter beings		
 DOM Series 		
+ ETCarren		
* GA Epitera		
* Mechanic Settings		
Michael Software Indian and Raccom		

Fig. 14: General Machine Settings UI

Machine Geometry Settings

The Geometry tab of the Machine Settings block allows the user to view and edit Machine (Equipment) coordinate system settings. These settings include the reference gantry angle, reference couch angle, reference collimator angle, and the rotation direction of each axis, relative to IEC 61217 coordinate system.

aninai Butonf		- I X
Mann Nas		
.dectred en		illeck Neper Kless (1941) 1240 🔒 🌅
 Paleni Saanh Saath la pálenk Irita spiere. 	Deneral Countries Australian Communities Advected	
 Import Patients 	Rati Reter. Octore Tech RetArge: 8-bg	
Y Department Configuration	Saudi Rublok. Disata Sudinas Calenais Rufingis: Educ	
+ Dynamics Relation	California: Reality: Counter Countries	
Countings		
Non Faculty Serlings Relation The Dirac additionations for a landing		
Mr. 1040-packed +		
• The belongs		
 DOM Settings 		
+ ETCarre		
• GA Epilona		
* Mechanic Settings		

Fig. 15: Machine Geometry Settings UI

7/9

sizes are displayed and the user can choose which applicators are available when creating a treatment

Users are unable to change manufacturer specific settings (such as physical block size).

Varian, or Elekta) chosen during the creation of a new machine.

The Commissioning tab of the Machine Settings block allows the user to view and edit commissioning energy data. The user can view the nominal energy and R90 values of the current commission data or import new commissioning data from a local file. The ElectronRT App currently only has support for commissioning data from the Pinnacle treatment planning system.

Commissioning Data Settings



The Applicators tab of the Machine Settings blocks allows for the viewing and/or editing of applicator and beamline device settings. Information regarding blocks, intensity modulators, and applicators of varying

plan using the current machine. The available applicators are dependent on type of machine (Siemens,

Applicator Settings

🔤 desireal Restored		- D X
Spin New		
Juctinal an		Allenia Nagar Khana 🛲 🔂 🛃
► Paleni Gash	Orenal Operative Systems Connections Advanced	
Saunt is palants in the system	Service Every 1 Sec	
* Import Patients	MR 187 (B	
Select DilliNipelants in Input	Minuting The	
¥ Ergeninelus Configuration	Kond Day	
+ Dynakalas kelings	RE 4DEas	
Exercise	His Little	
The later of the later	Kinish Dang 2000	
Natural action (1) to well wellings for a facelity	No SEVIN	
No: 1000adamit +	Inspected BL conversion data	
· to taken	Cellular Dec. 4311 SMD	
 DCBM Settings 	Output hacks make Application Taxa solution	
+ CTGames	Feerlike per een Tor's Tor's Dich State	
 GA Epilicia 		
	64 030106108100108	
* Machine Sellings	10-10 1001100110011001	
No. Ann Indings	2600 180 180 180 188	
tie er eit fastion	2020 100 100 100 100 100	
Veranlandmen		
	TWO BOA SHALTO THE THE THE THE	
	30 100 000 100 100 100	
	64 130 CM 1 CM 150 150	,
	10-10 1301 CHU CHU COU 1300	,
	Teris 1301 ENC/ ENC/ ESC 2000	
	2020 100100010000000	
	2624 1 801 6 801 6 801 6 801	
	E Hanina Dange O Mar	
	a spond (seg) to be	
	 Normal Loop: PLBAY 	
	Capit Leight M. Convention Table Tail	
	Destellinsern	

Fig. 17: Commissioning Data Settings UI

Advanced Settings

The Advanced tab of the Machine Settings block allows for the view and editing of advanced machine configuration parameters. These settings include tolerance tables, wedge tray information, blocking tray information, and room imaging parameters.

animal Redeved?		- п к
.dectred an		Allech Hayar Flores CHAR 1298: 😝 📳
P Andre Grenk Policie	Otherentin Space of the state	
 Exclusion ECHI Satispi ECHI Satispi ECHI Satispi ECHI Satispi ECHI Satispi 	Swalshink K BOX 7 BOX	-
A Martine Malayee		

Fig. 18: Advanced Machine Settings UI

