

# Prescription

You can add multiple prescriptions of different target structures and all of them will be applied to the treatment plan.

In the Prescription block you can set the following fields:

- **Fractions** : The number of fractions over which the prescription dose will be delivered.
- **Prescriptions** : the total prescription dose value in Gy (maximum of 100 Gy) and the target structure for said dose.
  - Note that at this time prescription values do not have an exact specified meaning (e.g. mean, min) as the ElectronRT application only uses these values to scale the relative dose to absolute dose in the UI.
- **Plan Re-Normalization** : The value used to re-normalize the plan level dose.




Fig. 1: Prescription Example

From:  
<http://apps.dotdecimal.com/> - **decimal App Documentation**

Permanent link:  
<http://apps.dotdecimal.com/doku.php?id=electronrt:userguide:tutorials:prescription&rev=1620758884>

Last update: **2021/07/29 18:24**

