



make this page much more useful and pretty

Thinknode Errors

When encountering an error in the astroid Planning App, check this list of common and known issues and their known fixes. If the fix detailed does not resolve the issue please contact the .decimal staff.

Error Messages

When receiving an error notification, this message will be in the body of the notification.

- **“provider terminated unexpectedly”:**
 - **During feasibility:**
 - Explanation: The calculation grid is too dense and crashed the provider.
 - Solution: Make calculation grid coarser by adding more ring structures and/or increasing the base resolution and region spacing

HTTP Error Codes

When receiving an error notification, this error code will be in the body of the notification.

- **401**
 - Explanation: The current Thinknode user session token has expired.
 - Solution: Close and reopen the Planning App, relogin to your realm, and re-initiate the task that failed.
- **409**
 - Explanation: The opened plan/patient was updated in the Thinknode record externally of the current Planning App instance. This resulted in the current instance being outdated.
 - Solution: Hit the F5 key or close and reopen the Planning App, relogin to your realm, and re-initiate the task that failed.
- **413**
 - Explanation: The calculation request body sent via the http request was too large to be accepted by Thinknode
 - Solution: Contact .decimal Staff. Please provide the complete error message and a full detailed description of the task attempted
- **50X**
 - Explanation: Any 500-504 error codes mean there was a temporary issue on the Thinknode server side.
 - Solution: Close and reopen the Planning App, relogin to your realm, and re-initiate the task that failed.

From:

<http://apps.dotdecimal.com/> - **decimal App Documentation**

Permanent link:

<http://apps.dotdecimal.com/doku.php?id=planning:userguide:errors&rev=1470765761>

Last update: **2021/07/29 18:22**

