



Launching Astroid

The Astroid *Launcher* will house the various applications that will be used as part of the Astroid Treatment Planning system. Any updates to these applications will automatically be deployed to the *Launcher*. The user will be notified that there is an update or a new version that will need installation once they have chosen an application.


1. Open the Astroid Launcher




2. Sign in using your user name and password

The image shows a "Sign In" form with a dark blue background. At the top, the text "Sign In" is centered in white. Below it, there are two input fields: "Username*" and "Password*", both with white text and a blue outline. The "Password*" field has a blue eye icon to its right. Below the input fields is a blue button with the text "LOGIN" in white.


3. Select the blue *Log In* button
4. Select your realm from the list of available realms

The image shows a "Realm Selection" form with a dark blue background. At the top, the text "Realm Selection" is centered in white. Below it is a list of realms, each with a name and a description. The realms are: "dev-anegron" (The development realm for Alfredo Negron-Rivera at .decimal.), "dev-aclark" (The development realm for Adam Clark at .decimal.), "dev-kmobile" (The development realm for Keitt Mobile at .decimal.), and "dev-staging" (The realm for the staging branch of work.). A red arrow points to the "dev-kmobile" realm.

5. From the left side tool bar select the application you would like to launch



6. Select the blue *Launch* button



7. If there is an updated version of this application the *Launch* button will say *Install*. Click this and allow the latest version to install. After the latest version is installed the blue button will say *Launch*
8. You may now proceed with opening a current patient or importing a new patient

From:

<http://apps.dotdecimal.com/> - decimal App Documentation

Permanent link:

<http://apps.dotdecimal.com/doku.php?id=planning:userguide:tutorials:launcher&rev=1484169835>

Last update: 2021/07/29 18:24

