


# Launching Astroid

The Astroid *Launcher* will house the various applications that will be used as part of the Astroid Treatment Planning system. Any updates to these applications will automatically be deployed to the *Launcher*. The user will be notified that there is an update or a new version that will need installation once they have chosen an application. Use the following steps to open the *Launcher* and launch the Astroid Planning App:

1. Open the *Launcher*




2. Sign in using your account, user name and password

The image shows a "Sign In" screen with a dark background. It has three input fields: "Account\*" with the value "decimal", "Username\*" with the value "kmobile", and "Password\*" with masked characters "\*\*\*\*\*". A blue "LOGIN" button is located at the bottom right.


3. Click the blue *Login* button
4. Select your realm from the list of available realms

The image shows a "Realm Selection" screen with a dark background. It displays a list of realms in a scrollable container. The realms are: "dev-anegron" (development realm for Alfredo Negron-Rivera at .decimal), "dev-aclark" (development realm for Adam Clark at .decimal), "dev-kmobile" (development realm for Keitt Mobile at .decimal), and "dev-staging" (the realm for the staging branch of work). A white arrow points to the "dev-kmobile" realm.

5. From the left side tool bar select the application you would like to launch



6. Click the blue *Launch* button



1. If there is an updated version of this application, the *Launch* button will not appear and instead an *Install* button will be available. Click this to allow the latest version to install. After the latest version is installed, the blue button will revert back to *Launch*, which you can now click

7. This will open the Astroid Planning application and bring you directly to the main patient search screen

8. You may now proceed with opening a current patient or importing a new patient

From:  
<http://apps.dotdecimal.com/> - **decimal App Documentation**

Permanent link:  
<http://apps.dotdecimal.com/doku.php?id=planning:userguide:tutorials:launcher&rev=1498574372>

Last update: **2021/07/29 18:24**

