


Plan History

The *Plan History* allows the user to see all changes and modifications made to a plan. Each modification to a plan will result in a new *Revision* of the plan. (note: in a complex plan this may result in hundreds of revisions) If the user is not satisfied with the results of a modification they may go back to any of the previous *Revisions* and make that instance the active instance. (note: this will cause another *Revision* to be created) Below is a walk through and explanation of each step.

1. Open the *General* block from the Plan Overview
2. Expand the *History* tree node. The number of plan revisions will be listed. In this below example there are 36 revisions
3. Select the *More Details* option in the bottom right hand corner



4. The summary of the latest revisions will expand to the right

| General | Revision 36: view make active |
|---------------------------------------|-------------------------------|
| 1 | created: 2016-10-04 09:21:34 |
| | by: Keitt Mobile |
| | summary: show |
| Clinical Guide | Revision 35: view make active |
| max of 76 Gy(RBE) to CTV | created: 2016-10-03 10:54:29 |
| min mean of 72 Gy(RBE) to CTV | by: Keitt Mobile |
| | summary: show |
| History | Revision 34: view make active |
| revisions: 36 | created: 2016-09-19 09:47:29 |
| initial revision: 2016-09-19 14:08:40 | by: Keitt Mobile |
| created by: Keitt Mobile | summary: show |
| latest revision: 2016-10-04 09:21:34 | Revision 33: view make active |
| modified by: Keitt Mobile | created: 2016-09-15 09:32:52 |
| Link Details | by: Keitt Mobile |
| | summary: show |
| Patient Geometry | Revision 32: view make active |
| 16 segmented structures | created: 2016-09-19 09:38:54 |
| 1 planning structure | by: Keitt Mobile |
| | summary: show |
| RSP Image | Revision 31: view make active |
| 0 structure(s) have density overrides | created: 2016-09-15 09:47:29 |
| | by: Keitt Mobile |
| | summary: show |
| Calculation Grid | Revision 30: view make active |
| 0 even Dose Spacing | created: 2016-09-15 10:07:22 |
| 0 Grid Regions | by: Keitt Mobile |
| | summary: show |
| Beams | Revision 29: view make active |
| 0 plan beams | created: 2016-09-16 08:44:49 |
| | by: Keitt Mobile |
| | summary: show |
| Optimization | Revision 28: view make active |
| 3 fraction groups | created: 2016-09-16 08:37:23 |
| 6 optimization constraints | by: Keitt Mobile |
| 1.4 optimization objectives | summary: show |
| | Revision 27: view make active |
| | created: 2016-09-15 09:36:56 |
| | by: Keitt Mobile |
| | summary: show |
| | Revision 26: view make active |
| | created: 2016-09-16 08:36:08 |
| | by: Keitt Mobile |
| | summary: show |
| | Revision 25: view make active |
| | created: 2016-09-15 09:21:09 |
| | by: Keitt Mobile |
| | summary: show |
| | Revision 24: view make active |
| | created: 2016-09-15 17:44:36 |
| | by: Keitt Mobile |
| | summary: show |
| | Revision 23: view make active |
| | created: 2016-09-15 17:44:12 |
| | by: Keitt Mobile |
| | summary: show |
| | Revision 22: view make active |
| | created: 2016-09-15 10:20:29 |
| | by: Keitt Mobile |
| | summary: show |

- The user should then pick the revision they want to look at by selecting the blue *show*. This will expand that particular version of the plan and explain what was done
 - “op”: is the operation performed in this version that differs from the previous version (i.e. delete: an item was deleted)
 - “path”: is the item in the treatment plan data that had the operation performed
 - “val”: is the value in the path that had the operation performed
- If the user would like to return to this version they may do so by selecting the blue *make active*
- The user can close the expanded version by selecting the blue *hide*

| |
|--|
| Revision 36: view make active |
| created: 2016-10-04 09:21:34 |
| by: Keitt Mobile |
| summary: hide |
| { |
| op: update, |
| path: [captured scope, directive], |
| val: { |
| some: 57f2736f00c0b9f2c378152e06cfd09 |
| } |
| { |
| op: update, |
| path: [captured scope, site], |
| val: { |
| some: 57dbebec00c0102cd41ea86de124b527 |
| } |
| } |

From:

<http://apps.dotdecimal.com/> - **decimal App Documentation**

Permanent link:

http://apps.dotdecimal.com/doku.php?id=planning:userguide:tutorials:plan_history&rev=1476737252



Last update: **2021/07/29 18:25**