

Plan Templates


Plan Templates are a powerful feature of Astroid that can greatly reduce the time spent in creating high quality treatment plans. In a simple sense, a *Plan Template* is just a starting point that can be used when creating new treatment plans for this or other patients. When a *Plan Template* is applied during new plan creation (see [here](#) for creating a new plan), all user specified plan parameters from the template are copied into the new plan. In essence you are applying the plan to a new set of images and structures, therefore it is important that structure names and other data are consistent in your *Patient Models* when using *Plan Templates*.

Saving Plan Templates

Plan Templates can be defined globally at the site level for all users at your institution and each user is also free to create their own private list of user level *Plan Templates*. A *Plan Template* can be easily created from any existing plan by opening the General/History block. At the bottom of the block the user may then choose either *Save as User Plan Template* or *Save as Site Plan Template*

The screenshot shows the Astroid planning interface. The top navigation bar includes the 'astroid' logo and the path 'planning > TG244 Prostate_Bed [002441] > Four Field'. Below this is a dropdown menu with 'General / History' selected. The main content area is divided into several sections: 'Plan Info' (plan name: Four Field, plan description: [empty], Edit / Rename, Publish), 'Clinical Goals' (no clinical goals), 'Prescriptions' (Implementing all course prescriptions, Edit Prescription), 'History' (revisions: 2, initial revision: 2018-09-25 10:57:26, created by: Keitt Mobile, latest revision: 2018-09-26 13:26:33, modified by: Keitt Mobile), 'IDs' (RKS: 5baa4cd601803baa58c5f86e91ce08e6 copy, ISS: 5babc14900c0c1ccb2dc3e63f8e8e2d6 copy), and 'Plan Templates'. At the bottom of the 'Plan Templates' section, there are two buttons: 'Save as User Plan Template' and 'Save as Site Plan Template'. A red arrow points from the 'Plan Templates' section towards these buttons.

The user then label and fill in any description information as needed.

A screenshot of a web application window titled "Plan Templates". Inside, there is a section titled "User Template Details". It contains two input fields: "Name:" with the text "Four Field" and "Description:" with a large empty text area. At the bottom right of the form are two buttons: "Save Plan Template" and "Cancel". Below the form, there is a link that says "Save as Site Plan Template".

When saving an existing plan as a template, there are a few rules/tips that should be kept in mind which help make your templates easier to use in the future-

- In order for ones templates to work well there needs to be consistency in ones Site Structure Templates, i.e.- Brainstem consistently labeled Brainstem, not brainstem, Brain stem, brain_stem etc.
- When using explicit points in templates be aware that they may not transfer as desired from plan to plan. Please check each explicit point on each plan
- Templates must be validated prior to the creating a new treatment plan to ensure compatibility with the new patient model
- If a template is created for a prostate patient that template will not show up in a breast patient model due to the patient structure model being different
- As a user one can create a template to any level into the planning process they would like- for example if one would like to stop at putting the calculation grid on they may or they may go as far as the optimization process

Using Plan Templates

Plan Templates are applied when creating a new plan:

1. Click on the blue *Add Plan* link under the *Patient Model* to create a new plan
2. From the Base Plan drop down menu choose the *Plan Template* that should be used and then click blue OK

The screenshot shows the 'astroid planning' application interface. At the top, there's a header with 'astroid' and 'planning > TG244 Prostate_Bed [002441]'. Below this, patient information is displayed: Name: TG244 Prostate_Bed, MRN: 002441, Sex: Male, Date of Birth: 1955-09-11. A checkbox for 'Show Archived Entries' is present. The main area shows a plan with Description: none and Prescription: Implementing all course prescriptions. A 'Plan:' dropdown is open, showing a list of templates: 'none (new plan)', 'Site Templates', 'User Templates', 'Current Patient', and 'Base Plan: none (new plan)'. The 'New Plan' section is highlighted, showing 'Name: opposed lats', 'Description: four field', and 'Base Plan: none (new plan)'. There are 'OK' and 'Cancel' buttons at the bottom right.

3. Open the new plan and check that all plan parameters have been applied as expected
4. Make any changes necessary within plan and then proceed with planning from this starting point
5. If used well, you should be able to review the plan, make a few minor edits, and be at the point of running feasibility & optimization within just a few minutes

From:
<http://apps.dotdecimal.com/> - **decimal App Documentation**

Permanent link:
http://apps.dotdecimal.com/doku.php?id=planning:userguide:tutorials:plan_templates&rev=1537986448

Last update: **2021/07/29 18:25**

