

Defining Constraints

About Constraints


Constraints can be set in multiple levels (*Fraction Group*, *Beam Set Group*, *Plan* and they will function differently in each level. *Constraints* applied at the *Fraction Group* level apply to all the beams as a whole regardless of the number of *Beam Sets*. Constraints applied at the *Beam Set* level apply to that group of beams. The *Constraint* dose is divided by the beams in the *Beam Set*. The following will provide a walk through of the different levels and how constraints work at each one. *Constraints* applied at the *Plan* level are considered “hard constraints”- constraints that have to be achieved.

Working with Constraints

Working with Fraction Group Constraints

The *Constraints* applied at the *Fraction Group* level take into account everything in the *Fraction Group* and *Beam Sets*

1. Select the *Fraction Group* if it has been created or create a new *Fraction Group*




2. Choose from the drop down the structure or structures that should be considered in the *Fraction Group*



3. Define what constraint(s) should be applied to each structure by choosing the constraint and entering the Gy



- Min- The minimum dose the structure should receive
 - Max- The maximum dose the structure should receive
 - Min Mean-the minimum mean dose a structure should receive
 - Max Mean- the maximum mean dose a structure should receive
 - the user can choose to apply one of these constraints or multiple constraints to the structure
4. Once the constraints are set in the *Fraction Group* the user will assign constraints in the *Beam Sets* associated to the *Fraction Group*. The assigned dose will be divided between the beams in the *Beam Set*



- Each individual *Beam Set* group can have constraints.
- The constraints only apply to the *Beam Set* group that they are associated with
- As before the user can apply more than one *Constraint* to a *Beam Set*
- Min- The minimum dose the structure should receive
- Max- The maximum dose the structure should receive
- Min Mean-the minimum mean dose a structure should receive
- Max Mean- the maximum mean dose a structure should receive

Working with Plan Constraints

The *Constraints* applied at the *Plan* level apply to the plan as a whole. These constraints are what drive the feasibility calculation- whether the plan is achievable. The calculation engine considers these Constraints as “hard constraints” that have to be achieved

1. Open the Constraints block contained in the Plan Constraints/Objectives block
2. Choose the structures to apply constraints to

3. Choose the constraints to apply to the structure and the Gy that needs to be achieved

Plan Constraints / Objectives

Constraints

PTV

Min: 36.25 Gy(RBE)

Max: 46 Gy(RBE)

Min Mean: add statement

Max Mean: add statement

Rectum

Max: add statement

Max Mean: 20 Gy(RBE)

s1 - Prostate + 0mm (in 3D)

Min: 40 Gy(RBE)

Max: 115 Gy(RBE)

Min Mean: add statement

Max Mean: add statement

Add Structure

OK Cancel

- Min- The minimum dose the structure will receive
- Max- The maximum dose the structure will receive
- Min Mean-the minimum mean dose a structure will receive
 - this will drive the dose up across the structure
- Max Mean- the maximum mean dose a structure will receive
 - this will limit the mean dose across the structure

Once all the Constraints have been set the user can either start the Feasibility or move on to defining the Objectives

From:
<https://apps.dotdecimal.com/> - **decimal App Documentation**

Permanent link:
https://apps.dotdecimal.com/doku.php?id=planning:userguide:walkthroughs:defining_constraints&rev=1470418432

Last update: **2021/07/29 18:25**