

# decimal Launcher Version History



This section is intended to provide access to previous release notes for each version of the decimal Launcher. Below you will find the userguide for each released version of the app. To find your application's current version use the top "Help" menu and select "About".

---

## Release Notes

### decimal Launcher 1.2.0

April 7, 2025

#### *What's New*

- Improved support for "All User" installations
- Updates can now be phased in and are not immediately forced
- Security updates

#### *Bug Fixes*

- Fixed hidden login information bug
  - Fixed error with Updater not installing over prior version
- 

### decimal Launcher 1.1.3

April 10, 2024

#### *What's New*

- Updated Auth token storage to fix 401 errors

#### *Bug Fixes*

- None
-

## decimal Launcher 1.1.2

March 10, 2023

### *What's New*

- None

### *Bug Fixes*

- HOTFIX: Revoked the JWT Token Check to resolve ssl certificate errors
- 

## decimal Launcher 1.1.1

February 13, 2023

### *What's New*

- Added a new check for JWT token validity before testing connection to the Direct API
- Added a right click Menu to the Apps list with "Open in explorer" and "Manage in Direct" (Manager only) options.
- Reworked "app is outdated" banners to be specific to managers if applicable, and made more clear.
- Added a refresh option under the app menu to manually refresh the Launcher application.
- Added a "Launching" state when an app is launched for clarity.
- Added a second alert by the Launch/ Download button that there is a new version of an app available for visibility.

### *Bug Fixes*

- None
- 

## decimal Launcher 1.1.0

October 03, 2022

### *What's New*

- Added alert banners to the Login page
- Reworked the error page to add better instructions and make the error more clear

### *Bug Fixes*

- None
- 

## decimal Launcher 1.0.6

May 25, 2022

### *What's New*

- Added New settings menu for changing the App timeout
- Sandbox mode now ignores apps that don't support a sandbox mode
- Added a better message for users that do not have any sites

### *Bug Fixes*

- Fixed a bug with Invalid app versions showing up as both invalid and out of date
- 

## decimal Launcher 1.0.5

December 15, 2021

### *What's New*

- Adds a refresh button to the application list to manually refresh apps and versions after they've been changed in decimal Direct

### *Bug Fixes*

- Fixes a bug that didn't require MFA for specific users when MFA is forced
- 

## decimal Launcher 1.0.4

September 23, 2021

### *What's New*

- None

### *Bug Fixes*

- Fixes a bug where MFA would not save when using the 'Remember this Device' option
  - Fixes a bug with sites not showing up in the correct order in the site list
-

## decimal Launcher 1.0.3

July 27, 2021

### *What's New*

- Adds downloads from Azure instead of .decimal servers to allow for faster downloads
- Adds an icon to link to the user preferences page to decimal Direct

### *Bug Fixes*

- Fixes a bug that would result in an endless loop when attempting to download an app
- 

## decimal Launcher 1.0.2

June 17, 2021

### *What's New*

- Adds site nicknames to the site ID list
- Adds a progress bar when downloading an application

### *Bug Fixes*

- Fixes application downloads failures putting the application in a bad state
- 

## decimal Launcher 1.0.1

January 12, 2021

### *What's New*

- Adds a link for new application versions to automatically send the user to decimal Direct's manage page
- Adds a manual download option for if the automatic download fails

### *Bug Fixes*

- Fixes an error checking for a connection to .decimal servers when logging in
  - Fixes a bug with applications that have a space in their executable name
-

# decimal Launcher 1.0.0

August 19, 2020

## *What's New*

- Initial Release of the decimal Launcher

## *Bug Fixes*

- None
- 

From:

<http://apps.dotdecimal.com/> - **decimal App Documentation**

Permanent link:

<http://apps.dotdecimal.com/doku.php?id=decimallauncher:versions&rev=1744034983>

Last update: **2025/04/07 14:09**

