

# Launching ElectronRT

In order to Launch the decimal ElectronRT app you must use the decimal Launcher to authenticate and open your desired application.

**NOTE:** You will need your NEW Auth0 Direct user Log in.

- Download and install the decimal Launcher.
- Log in with your Auth0 credentials.
- Find decimal ElectronRT and download/launch the application.

For further instructions please refer to the decimal Launcher User Guide. [decimal Launcher Overview](#) > User Guide

From:  
<http://apps.dotdecimal.com/> - decimal App Documentation

Permanent link:  
[http://apps.dotdecimal.com/doku.php?id=electronrt:userguide:tutorials:launching\\_electronrt&rev=1594561514](http://apps.dotdecimal.com/doku.php?id=electronrt:userguide:tutorials:launching_electronrt&rev=1594561514)

Last update: 2021/07/29 18:24

