

Launching ElectronRT

In order to Launch the decimal ElectronRT app you must use the decimal Launcher to authenticate and open your desired application.

NOTE: You will need your NEW Auth0 Direct user Log in.

1. Download and install the decimal Launcher.
2. Log in with your Auth0 credentials.
3. Find decimal ElectronRT and download/launch the application.

For further instructions please refer to the decimal Launcher User Guide. [decimal Launcher Overview](#) > User Guide

From:
<http://apps.dotdecimal.com/> - **decimal App Documentation**

Permanent link:
http://apps.dotdecimal.com/doku.php?id=electronrt:userguide:tutorials:launching_electronrt&rev=1594561852

Last update: **2021/07/29 18:24**

