

p.d Version 6 (with OptiBol)



Version 6 is the first major software release of p.d in over a decade. Version 6 maintains the same look and feel as the previous version 5 software, but is greatly enhanced by the three major updates described below.

OptiBol™

Built on the Benefits of FlexiBol®, Optimized for Today's Radiotherapy Workflows

- Automated, variable-thickness bolus saves you time in bolus design and treatment planning
- Target-based design ensures the dosimetric buildup you need, where you need it, for any potential photon beam geometries
- Built-in alignment guides serve as an option to help radiation therapists in bolus and patient setup
- No added cost means the same device pricing as FlexiBol, with access for all customers

[EXPLORE OPTIBOL ON YOUR OWN](#)

[REQUEST A PERSONALIZED DEMO](#)

Bolus Wizard (Guided Approach)

- Major overhaul to Bolus Wizard workflow to improve clarity and enhance the end user experience
- Adds Bolus Size and Price information
- Adds checks for common errors/inconsistencies (e.g. check for unintended holes and disconnected pieces)
- Adds support for ordering Brass Electron Skin Collimators (please call to request access to this option)

Enhanced Order Wizard

- Expected Ship Date and Estimated Delivery Date/Time now shown when selecting order priority

From:

<http://apps.dotdecimal.com/> - **decimal App Documentation**

Permanent link:

http://apps.dotdecimal.com/doku.php?id=pdotd:pdotd_version_6&rev=1782752424

Last update: **2026/06/29 17:00**

