

# p.d Version History

This section is intended to provide access to previous versions of the p.d documentation for releases of p.d after 5.4.0.

---

## Release Notes

### p.d 5.4.0

July XX, 2021 -

#### *What's New*

- Updated p.d to be installed and launched from the decimal Launcher application
  - Removed the login screen as p.d must always be launched using the decimal Launcher
  - Removed internal update mechanics from p.d as install and updates will be managed by the decimal Launcher
- Updated the Import Wizard machine drop downs to only include machines with the current application site ID (sent from the decimal Launcher)
- Added a new Planning Menu option to take the user to the user profile management page of decimal Direct

#### *Bug Fixes*

- none
- 

From:  
<http://apps.dotdecimal.com/> - **decimal App Documentation**

Permanent link:  
[http://apps.dotdecimal.com/doku.php?id=pdotd:pdotd\\_versions&rev=1627304820](http://apps.dotdecimal.com/doku.php?id=pdotd:pdotd_versions&rev=1627304820)

Last update: **2021/07/29 18:19**

