

# p.d Version History



This section is intended to provide access to previous versions of the p.d documentation for releases of p.d after 5.4.0.

---

## Release Notes

### p.d 5.4.1

September 20, 2021

#### *What's New*

- Updated all design wizard machine drop downs to only include machines for the current site facility (selected from the decimal Launcher)
- Updated p.d to no longer include or uses addresses from the treatment machine
- Updated p.d ordering address to use the decimal Direct shipping address and disabled changing of the ordering address in the Order Wizard

#### *Bug Fixes*

- none

---

### p.d 5.4.0

August 26, 2021

#### *What's New*

- Updated p.d to be installed and launched from the decimal Launcher application
  - Removed the login screen as p.d must always be launched using the decimal Launcher
  - Removed internal update mechanics from p.d as install and updates will be managed by the decimal Launcher
- Updated the Import Wizard machine drop downs to only include machines with the current application site ID (sent from the decimal Launcher)

- Added a new Planning Menu option to take the user to the user profile management page of decimal Direct

### *Bug Fixes*

- none
- 

From:

<https://apps.dotdecimal.com/> - **decimal App Documentation**

Permanent link:

[https://apps.dotdecimal.com/doku.php?id=pdotd:pdotd\\_versions&rev=1632148379](https://apps.dotdecimal.com/doku.php?id=pdotd:pdotd_versions&rev=1632148379)

Last update: **2021/09/20 14:32**