

# Astroid Patient Data Model

The following page describes the hierarchy of data used to manage patient data records within the Astroid planning environment.

## Hierarchy

- [illegible]

- Beamset [1,...,N]
  - Constraint [0,1,...,N]
  - Beam [1,...,N]
- Constraints [0,1,...,N]
- Objectives [0,1,...,N]
- Dose Results

From:

<http://apps.dotdecimal.com/> - **decimal App Documentation**

Permanent link:

[http://apps.dotdecimal.com/doku.php?id=planning:instructions\\_for\\_use:instructions\\_for\\_use&rev=1459184100](http://apps.dotdecimal.com/doku.php?id=planning:instructions_for_use:instructions_for_use&rev=1459184100)

Last update: **2021/07/29 18:22**

