

- Target
- Constraint [0,1,...,N]
- Beamset [1,...,N]
 - Constraint [0,1,...,N]
 - Beam [1,...,N]
- Constraints [0,1,...,N]
- Objectives [0,1,...,N]
- Dose Results

Descriptions

- Patient
 - A person receiving medical treatment.
- Course:
 - A prescribed regimen to be followed for a specific period of time.
- Intent:
 - The purpose of the treatment.
- Directive:
 - The prescription for the course of treatment.
- Snapshot:
 - Overview of targets and organs at risk.
- Request:
 - The date of implementation of the number of fractions. Can be all fractions or a specific number of fractions.
- Plan:
 - All aspects of the patient planning take place here. Such as patient geometry, calculation grid, beams and optimization.

From:
<http://apps.dotdecimal.com/> - **decimal App Documentation**

Permanent link:
http://apps.dotdecimal.com/doku.php?id=planning:instructions_for_use:instructions_for_use&rev=1464019275

Last update: **2021/07/29 18:22**

