2025/12/14 21:20 1/1 Planning Version History

## **Planning Version History**

This page is intended to provide access to previous versions of the astroid userguide documentation. Below you will find the userguide for each released version of the astroid Planning App.

## Finding the App Version

To find your application's current version:

- 1. Open the astroid Application
- 2. Login with your username and password
- 3. Select your realm
- 4. From the top View menu select App Info
- 5. The "Version:" field will contain your current version of the App

Using the current version of the application, find the appropriate userguide from below

## **Versions**

• 1.0.0-a11

(Alpha release 11)

From:

http://apps.dotdecimal.com/ - decimal App Documentation

Permanent link:

http://apps.dotdecimal.com/doku.php?id=planning:planning\_versions&rev=1474554079

Last update: 2021/07/29 18:19

