

Astroid Planning Release Notes

Version: 2.1.2

Date: Mar 29, 2018

The release version build of the Astroid Planning App version 2.1.2 is now available (note: this version is still pending FDA 510k clearance). No data types or api function signatures are changed in this release. This release includes the following changes since the 2.1.1 release:

- Fixes Results API bug with unsupported beam_isocenter union type

From:

<http://apps.dotdecimal.com/> - **decimal App Documentation**

Permanent link:

http://apps.dotdecimal.com/doku.php?id=planning:release_notes:2.1.2&rev=1572619836

Last update: **2021/07/29 18:22**

