

# Astroid Planning Release Notes

**Version:** 2.3.3

**Date:** Jun 3, 2019

The release version build of the Astroid Planning App version 2.3.3 is now available. No existing ResultsAPI function signatures are changed in this release. This release includes the following changes since the previous release:

- Added an option to apply RSP overrides to the CT image in the DICOM Export Block
- The ResultsAPI Viewer DICOM AE export servers now pull from the live site info in the realm

From:

<http://apps.dotdecimal.com/> - **decimal App Documentation**

Permanent link:

[http://apps.dotdecimal.com/doku.php?id=planning:release\\_notes:2.3.3](http://apps.dotdecimal.com/doku.php?id=planning:release_notes:2.3.3)

Last update: **2021/07/29 18:23**

