Astroid Planning Release Notes

Version: 3.2.0 **Date:** Feb 15, 2021

The release build of the Astroid Planning App version 3.2.0 is now available. No existing ResultsAPI function signatures are changed in this release, however the output has changed for functions returning dose or fluences due to the addition of the new optional fluence overrides. This release includes the following changes since the previous release:

- 1. Adds an optional<fluence override> to the treatment plan data ("overridden fluence")
 - 1. If present, all dose results will utilize these fluences instead of those computed by the Astroid Planning App
 - 2. Note, for the MGH Build, Low MU Filtering will be applied to the override fluences
- 2. (MGH Build Only) Applies Low MU Spot Filtering function to fluence results when not actively navigating the MCO results

From:

http://apps.dotdecimal.com/ - decimal App Documentation

Permanent link:

http://apps.dotdecimal.com/doku.php?id=planning:release_notes:3.2.0&rev=1613161701

Last update: 2021/07/29 18:22

