



make this page much more useful and pretty

# Thinknode Errors

- **“provider terminated unexpectedly”:**
  - During feasibility: calculation grid is too dense and crashed the provider. Make calc grid coarser
- **“401”:**
  - Session token has expired, close planning and relogin.
- **“50X”**
  - Any 500-504 error codes mean there was a temporary issue on the Thinknode side. Closing, reopening, and retrying what gave the 50X error should resolve the issue.

From:

<http://apps.dotdecimal.com/> - **decimal App Documentation**

Permanent link:

<http://apps.dotdecimal.com/doku.php?id=planning:userguide:errors&rev=1470765167>

Last update: **2021/07/29 18:22**

