

Plan History

The *Plan History* allows the user to see all changes and modifications made to a plan. Each modification to a plan will result in a new *Revision* of the plan. (note: in a complex plan this may result in hundreds of revisions) If the user is not satisfied with the results of a modification they may go back to any of the previous *Revisions* and make that instance the active instance. (note: this will cause another *Revision* to be created.

From:

<http://apps.dotdecimal.com/> - **decimal App Documentation**

Permanent link:

http://apps.dotdecimal.com/doku.php?id=planning:userguide:tutorials:plan_history&rev=1476714436

Last update: **2021/07/29 18:25**

