2025/11/21 15:34 1/1 Overview

## **Overview**

The .decimal astroid app suite is a system of applications living on the thinknode $^{\text{m}}$  framework. Each app has it's own intended use, documentation, and has been thoroughly tested before being released into the production environment.

Access to each App is provided by the thinknode<sup>TM</sup> framework using http json formatted requests. thinknode<sup>TM</sup> provides the 'backbone' used to send and receive requests, maintain users, realms and organizations, and provide data storage and management.

## **Astroid Apps**

Each of the released apps contains its own use documentation. Select the app below to see its corresponding user documentation.

**Dosimetry App** Planning and analysis of proton radiation therapy treatments. Includes design task, dose calculation, and radiotherapy support functions.

From:

http://apps.dotdecimal.com/ - decimal App Documentation

Permanent link:

http://apps.dotdecimal.com/doku.php?id=start&rev=1435851513

Last update: 2021/07/29 18:18

