

# Overview

The .decimal astroid app suite is a system of applications living on the thinknode™ framework. Each app has its own intended use, documentation, and has been thoroughly tested before being released into the production environment.

Access to each App is provided by the thinknode™ framework using http json formatted requests. thinknode™ provides the 'backbone' used to send and receive requests, maintain users, realms and organizations, and provide data storage and management.

## Astroid Apps

Each of the released apps contains its own use documentation. Select the app below to see its corresponding user documentation.



**Dosimetry App** Planning and analysis of proton radiation therapy treatments. Includes design task, dose calculation, and radiotherapy support functions.



**Dicom App** Parsing and storing of DICOM file types into standard data types. Includes support for Plan, Structure Set, CT Image, and Dose dicom files.

From:

<http://apps.dotdecimal.com/> - **decimal App Documentation**

Permanent link:

<http://apps.dotdecimal.com/doku.php?id=start&rev=1436370063>

Last update: **2021/07/29 18:18**

