2025/11/21 15:38 1/1 Overview

Overview


The .decimal astroid app suite is a system of applications living on the thinknode™ framework. Each app has it's own intended use, documentation, and has been thoroughly tested before being released into the production environment.

Access to each App is provided by the thinknodeTM framework using http json formatted requests. thinknodeTM provides the 'backbone' used to send and receive requests, maintain users, realms and organizations, and provide data storage and management.

Astroid Apps

Each of the released apps contains its own use documentation. Select the app below to see its corresponding user documentation.

Dosimetry App Planning and analysis of proton radiation therapy treatments. Includes design task, dose calculation, and radiotherapy support functions.



From:

http://apps.dotdecimal.com/ - decimal App Documentation

Permanent link:

http://apps.dotdecimal.com/doku.php?id=start&rev=1436370063

Last update: 2021/07/29 18:18

