2025/11/21 15:38 1/2 Overview

Overview


The .decimal astroid app suite is a system of applications living on the thinknode™ framework. Each app has it's own intended use, documentation, and has been thoroughly tested before being released into the production environment.

Access to each App is provided by the thinknode^{\dagger} framework using http json formatted requests. thinknode^{\dagger} provides the 'backbone' used to send and receive requests, maintain users, realms and organizations, and provide data storage and management.

Astroid Apps

Each of the released apps contains its own use documentation. Select the app below to see its corresponding user documentation.

Dosimetry App Planning and analysis of proton radiation therapy treatments. Includes design task, dose calculation, and radiotherapy support functions.



Manifest Documentation

The Manifest Guide for each app contains a complete list of all functions and types supported by the app through thinknode^{$^{\text{TM}}$} API calls. Refer to each astroid Apps Examples section for outlined usage on function calling and calculation requests as well as connecting to the thinknode^{$^{\text{TM}}$} framework.

Refer to the thinknode[™] links on the left sidebar for detailed help using the thinknode[™] framework.

Support

For questions, comments, or to schedule a training session, please contact our customer support team at: appsupport@dotdecimal.com

2025/11/21 15:38 2/2 Overview

Copyright © 2015 .decimal, LLC. 121 Central Park Place, Sanford, FL 32771. All Rights Reserved. astroid $^{\text{\tiny TM}}$ is trademark of .decimal, LLC. thinknode $^{\text{\tiny TM}}$ is trademark of Thinknode Labs, LLC.

From:

http://apps.dotdecimal.com/ - decimal App Documentation

Permanent link:

http://apps.dotdecimal.com/doku.php?id=start&rev=1439825538

Last update: 2021/07/29 18:18

