

# .decimal Legacy Account Migration

This article provides a detailed comparison overview of the .decimal Legacy account features and the new .decimal Auth0 account features. It also outlines, at a high level, the process which users will take to migrate their existing Legacy accounts to the new .decimal account.

## Account Type Comparison

The below table details the main difference in features between the Legacy .decimal accounts and the new .decimal Accounts.

| Feature  | Legacy                                     | Auth0 |
|--|--|-------|
| <b>Self Management</b>   |  |       |
| Self Change Password Option  | Yes  | Yes   |
| Forgot Password / Reset  | [.d]                                       | Yes   |
| Enable <a href="#">MFA/2FA</a>   |  | Yes   |
| Forced non-shared logins (unique to each user)                                   |  | Yes   |
| <b>Manager Level Account Management</b>  |  |       |
| Add & Invite Other Users   | [.d]                                       | Yes   |
| Remove users from site   | [.d]                                       | Yes   |
| Require all users to use <a href="#">MFA</a>                                     |  | Yes   |
| Push app version updates to all users  |  | Yes   |
| Change default shipping address in decimal Direct                                |  | Yes   |
| <b>Application Features</b>  |  |       |
| Use <a href="#">decimal Launcher</a> to install and launch .decimal applications |  | Yes   |
| Automatic email alerts for new application versions and releases                 |  | Yes   |
| <b>Legend Key</b>  | <b>Description</b>                         |       |
| [.d]   | .decimal staff must perform this task      |       |
| Yes  | This feature is available to end users     |       |
|  | This feature is not available to end users |       |

## Account Migration Overview

Currently .decimal staff are required to migrate your existing Legacy accounts to the new .decimal account. Once your existing accounts have been migrated your Site Manager will be able to [manage your site](#) or [manage your users](#).

## Migration Process

At a high level, the migration process for existing sites and users looks as follows:

## Account Migration

1. .decimal staff will send the new system and security requirements documents to you and your clinic (e.g.: p.d 5.4's [RN-29](#) and [RN-30](#))
2. Your clinic's IT department should ensure the new network requirements are satisfied
3. .decimal staff will prepopulate for your account:
  1. Your [Site Information](#)
  2. Add the licensed .decimal applications to your site
4. .decimal staff will conduct a quick site survey to determine who, if anyone, at your site should be given [Site Management](#) permissions
5. .decimal staff will migrate any desired existing Legacy accounts that should be converted.
  - Existing Legacy accounts will have their login username changed to their email address associated with the Legacy account. **Note: Users will be unable to use the Legacy login once this step has been performed.**
  - During the migration step each user will receive an email notification to the email address associated with their existing Legacy account
  - The email notification will contain a link for the existing Legacy user to set their password to complete the migration

## New Software Installation

1. Legacy users will be instructed to manually uninstall p.d versions from their local work stations
2. Each user at each work station will be required to login to [decimal Direct](#) and [download and install decimal Launcher](#)
  - This step should only take a few minutes once the user has already setup their new password from the previous step
3. Users can then open the decimal Launcher application, select their site (if your organization has more than one), and launch the appropriate .decimal application

From:  
<https://apps.dotdecimal.com/> - **decimal App Documentation**

Permanent link:  
[https://apps.dotdecimal.com/doku.php?id=support:legacy\\_account\\_migration&rev=1659644810](https://apps.dotdecimal.com/doku.php?id=support:legacy_account_migration&rev=1659644810)

Last update: **2022/08/04 20:26**